

MY75S (Malt & Yeast Extract - 75% Seawater Agar)

Marine Protozoa

Medium	per litre
Natural filtered seawater	750.0 ml
Deionised water	250.0 ml
Malt extract	0.1 g
Yeast extract	0.1 g
Bacteriological agar	15.0 g

Make up the 75% seawater and heat (do not bring to boil). Then, add and dissolve the malt extract and yeast extract. Cool slightly then add the agar and shake to disperse. Autoclave at 15 psi for 15 minutes.

Reviewed: 17th January 2025